



Labs21 Advanced Course Series

High-Performance Laboratory Exhaust Devices

Tom Smith, P.E., M.S.
Exposure Control Technologies

Geoffrey C. Bell, P.E., M. Arch.
Lawrence Berkeley National Laboratory

Goal: Select Advanced Exhaust Devices

Objectives: At the end of this session, you will be able to:

- **Identify different types of exhaust devices**
- **Describe the process to select appropriate devices for a given context**
- **Distinguish between standard, good and better practice in selecting devices and assuring their performance**

Outline

- **Overview**
- **Design Considerations**
- **Device Review**
- **Selection Process**
- **Design Practice**
- **Performance Examples**
- **Conclusion**

Overview

High-Performance System

- **SAFE-** Meets operating specifications and performance requirements.
- **DEPENDABLE** - Predictably performs to meet user demands
- **EFFICIENT-** Minimum energy use and lowest possible operating costs
- **FLEXIBLE** - Adaptable to changing research needs

Design Considerations

- **Primary Safety Device**

- Intended to protect operator(s) health and welfare
- Selected specifically for scientific procedures

- **HVAC System Interface**

- Influences size and life-time performance of a laboratory's heating, ventilation and air conditioning (HVAC) system
- Verify performance and proper installation during final test, balance, and commissioning

- **Airflow Requirements**

- Ensure device containment performance prior to occupying laboratory, especially in “high fume-hood-concentration” labs
- Verify adequate airflow for either internal load demand or air change rate

Design Considerations

Needs and Opportunities...

- **Improve Containment Capability**
- **Promote Safe Use and Ergonomics**
- **Reduce Required Flow / Energy Consumption**
- **Integrate Monitoring and Control for Dependable Operation**
- **Minimize Maintenance**

Design Considerations

- **Coping with hazardous pollutants**
 - Potential pollutants
 - Sources of pollutants
 - Spill scenarios; device application
 - User/operator skills and experience
 - Egress pathways; device location limitations
- **Meeting space-conditioning requirements**
 - Lab layout/arrangement and equipment
 - Supply diffuser type and placement
 - Return grill location and quantity
 - Differential pressure; airflow between rooms; isolation protection
 - Diversity: factor assessment
 - Hood exhaust flow versus general exhaust flow

Design Considerations

- **Codes**

- UBC, UMC, IBC, IMC
- CAL/OSHA 5154.1

- **Standards**

- ANSI/AIHA Z9.5 - 2003
- NFPA-45 - 2000
- ACGIH - Industrial Ventilation - 24th Ed. - 2001
- ASHRAE Laboratory Design Guide- 2001
- OSHA - 29 CFR - Part 1910.1450
- Scientific Equipment & Furniture Association, SEFA 1.2-2003

Design Considerations



Design Considerations



Know the operating limitations!

Device Review




Types of Laboratory Exhaust Devices...

- **Chemical Fume Hoods**
 - Bench-Top, Distillation, Floor mounted (a.k.a. Walk-in), Scale-up
 - Radioactive, Perchloric, Auxiliary Air
 - High Performance, VAV, etc.
- **Biological Safety Cabinets**
 - Class I, Class II, Class III
 - Type A, Type B1, Type B2
- **Laminar Flow Fume Hoods**
- **Balance/Weighing Hoods**
- **Canopy Hoods**
- **Snorkel Exhausts**
- **Glove Boxes**
- **Microelectronic Hoods**
- **Ventilated Enclosures**

Device Review

- **Types of Laboratory Hoods**

- Conventional Fume Hood
- Constant Volume Bypass Fume Hood
- Variable Air Volume Fume Hood
- Restricted Opening Low-Flow Fume Hood
- Low-Velocity Low-Flow (LV/LF) Fume Hood



Device Review

Vertical Sash Opening

- Most common sash
- Good horizontal access
- Best with sash stop

Vertical
Sash Stop




Device Review

Horizontal Sash Opening

- Can reduce airflow volume
- May increase worker safety
- Caution – sash panels can be removed; defeats safety

Sash Panels



Device Review

Auxiliary Air Hood...

- Not Recommended
- Reduces containment performance
- Decreases worker comfort
- Disrupts lab temperature and humidity




Device Review

Distillation Hood

- Specific use
- Convenient experiment setup
- Optimizes experiment observation

Dual Sash Panels



Device Review

Floor Mounted Hood (Walk-In Hood)

- Specific use
- Optimizes experiment observation
- Limited containment performance

Horizontal Sliding Sash Panels



Device Review

Ductless Fume Hood

- Particulate Use
- Limited for Gas & Vapors
- Not Recommended for Permanent Installation



Device Review

Weighing Enclosure – Balance Hood

- Low-hazard use
- Minimal airflow
- Small footprint
- Requires exhaust connection



Device Review

Large Canopy Hood

- Heat Removal
- Energy Hog



Device Review

Exhaust Snorkel Hoods

- Low-hazard application
- Tasked device for vapors, e.g. soldering
- Close damper when not in use



Device Review

Ventilated Enclosure


- Minimal airflow volume
- Infrequently accessed
- Requires external exhaust
- Can reduce lab hazard classification



Device Review


- **Variable Air Volume Fume Hoods:**
 - Decrease laboratory exhaust flow rate
 - Decrease conditioned supply air
 - Reduces heating and cooling costs
 - Reduces supply and exhaust fan horsepower
 - Monitor sash position remotely
 - Check face velocity and laboratory pressurization
 - Provided constant face velocity control

Device Review




VAV Hood Operation

Device Review



VAV System

Device Review



VAV System

Device Review

- **Does the Low Flow / Low Velocity Hood provide:**
 - Equivalent or Better Containment at Reduced Face Velocities and Flow Volumes?
 - Improved performance for all users, even under misuse conditions?
 - More Robust and Less Susceptible to External Factors?
 - Better Monitoring and Flow Control?

If so... = High Performance Hood

Device Review

- **High-Performance Hood:
Improved Performance Through Better Design...**
 - Aerodynamic Entry
 - Directed Air Supply
 - Perforated or Slotted Rear Baffle
 - Airfoil Sill and Sash Handle
 - Integrated Monitors
 - Interior Dimensions

Device Review

- **High Performance Fume Hoods: current fabricators...**
 - Lab Crafters
 - Labconco
 - Fisher Hamilton
 - Kewaunee Scientific
 - Laboratory Equipment Manufacturers
 - Berkeley Hood
 - Others

Device Review


Lab Crafters Air Sentry HPFH

Upper chamber Turning Vane ●

Aerodynamic Sash Frame ●

Side Post Airfoils ●

Multi-Slot Front Airfoil ●



Device Review

Labconco XStream Hood

Modified Aerodynamic
Sash Pull

Modified Baffle
and Slots

Aerodynamic
Airfoil



Device Review

Fisher Hamilton PIONEER

- Automatic sash closer
- Directed supply flow @ full open sash
- Flush Airfoil Sill



Device Review

Berkeley Hood by LBNL

- Air Divider Technique
- Perimeter Air Supply
- Perforated Rear Baffle
- Slot Exhaust
- Optimized Upper Chamber
- Designed to minimize escape by reducing reverse flow



Device Review


Zone Occupancy Sensor

Sash Sensor/Monitor



Selection Process

Basic approach for safe and efficient results...



Selection Process...

- **Understand Hazards and Exhaust Device Use**
 - Type of Hazards
 - Frequency of Use: Dose
 - Duration of Use: Dose
 - Hazard Generation Characteristics
 - Effluent Characteristics

Selection Process...

- **Select Appropriate Laboratory Exhaust Device**

- Type
- Size
- Quantity
- Manufacturer & Model
- Operational Specifications
- Constraints and Limitations

**Re-check Hood and System is
suitable for the intended application**

Selection Process...

- **Key Factors Affecting Laboratory Exhaust Device Performance...**
 - Design features and installation
 - Lab design and layout
 - Ventilation system design (diffusers) and operation
 - Worker population and practices
 - Traffic patterns
 - Proximity to room egress
 - Laboratory equipment operation

Standard Design Practice

- **Choose Device...**

- Number
- Placement
- Features (verify with user)
- Chemical Resistance

- **Commissioning Requirements...**

- Volume Flow Check (TAB)
- Face Velocity Testing (per authority having jurisdiction)
- Certification (if required)
- EMCS tie-in (if available)

Good Design Practice


- **Choose efficient device...**
 - VAV hood Low-Flow hood
 - Consider alternative devices
- **Consider aspects of air management...**
 - Supply air temperature rate of change
 - Diffuser type and air “throw”
- **Require full ANSI/ASHRAE 110 testing as installed...**
 - Tracer gas containment
 - Sash movement effect
- **Conduct Operational Assessment Commissioning...**
 - Room differential pressure tests
 - VAV hood recovery tests (from unusual user interface, e.g., door slams)

Better Design Practice

- **Perform Computational Fluid Dynamic (CFD) Modeling...**
 - Model hood containment
 - Simulate lab space during spill event
- **Conduct advanced containment tests...**
 - Human-as-Mannequin (HAM) tests
 - Walk-up, walk-by tests
- **Require Performance Measurement Commissioning...**
 - Simulate “stress conditions”, e.g., whole building power failure
 - Test entire facility’s “harmony”
- **Arrange Efficiency Assurance Commissioning...**
 - Verify energy efficiency measures
 - Seasonal, full occupancy check


Performance Example: Good Practice

Observing reverse airflow
during ANSI/ASHRAE smoke
tests




Performance Example: Good Practice

ASHRAE 110 Tracer Gas Test



Center position



Right position

Performance Example: Better Practice

Preparing for Human-as-Mannequin tests

Tracer gas
detector
mounted
on human



Objects to be
manipulated




Performance Example: Better Practice

Conducting Human-as-Mannequin tests




Performance Example: Better Practice



CFD Fume Hood Models



Two-dimensional CFD



Three-dimensional CFD



Performance Example: Better Practice

CFD Lab Models

Princeton University

Conclusion

- **Design Mission a Success...?**
 - Increased operational safety significantly
 - Advanced Ergonomics
 - Minimized effects from external factors
 - Provided remote Monitoring and Control
- **Performance bottom line...**
 - Improved Containment
 - At Face Velocities of 50-60 fpm
 - With Full Range of Sash Opening
 - Independent of User Height

Conclusion

- **Summary Recap:**

- Design Considerations
- Device Review
- Selection Process
- Design Practice

- **Primary Issues:**

- Safety – Crucial reason for lab exhaust devices – Removal of Hazard.
- Temperature and humidity control – heat gain from equipment, computers, people.
- Productivity of facility – support mission.
- Cost to Design; to Build; to Operate.

Conclusion

It's a Combined Effort – A sustainable, energy-efficient facility requires teamwork involving:

- A “user” representative that understands the scientific research.
- A facilities representative that understands the building's HVAC systems.
- An industrial hygiene representative that understands worker health and safety issues.
- A laboratory design representative that understands exhaust device designs and features.

For More Information

Main Labs21 web site:

<http://www.labs21century.gov>

Primary Contacts:

Geoffrey C. Bell, P.E.
Lawrence Berkeley National Laboratory
Phone: 510 486-4626
Fax: 510 486-4089
E-mail: gcbell@lbl.gov

Otto Van Geet, P.E.
National Renewable Energy Laboratory
Phone: 303 384-7369
Fax: 303 384-7330
E-mail: otto_vangeet@nrel.gov